

# Vacuum Oven

## TVO-2-SV

PART ID: CTVS222

# CASCADE™ TEK

Capacity	1.67 cu.ft. (47.2 liters)
Interior Dimensions (W x D x H)	12" x 20" x 12" 304 x 508 x 304 mm
Exterior Dimensions (W x D x H)	24.5" x 34.9" x 27.1" 623 x 886 x 689 mm
Unit Weight	183 lb (83 kg)
Shipping Weight	224 lb (102 kg)
Shelves	3 standard, 3 total
Maximum Weight Per Shelf	35 lb / 15.8 kg
Permitted Total Load	105 lb / 47.6 kg
Standard Door Gasket	Silicone
Access Port	1" vacuum port (KF-25)
Controls	Watlow EZ-ZONE
Warranty	2 Year Limited Warranty
Interior Construction	300 series stainless steel
Exterior Construction	Painted Cold Roll Steel
Wall Clearance, Each Side	6.0" (152 mm)
Wall Clearance, Rear	12.0" (305 mm)
Vertical Clearance	12.0" (305 mm)
<b>Electrical, 50/60Hz, AC, 1Ø</b>	
Voltage	110V-120V
Full Load Amps	10
Nominal Power (watts @150°C)	662
Plug Supplied	NEMA 5-15P
<b>Energy Consumption (kWh/day)</b>	
At 80°C	12.8
At 150°C	15.4
At 220°C	24.4
<b>Temperature Range 10°C above ambient to 220°C</b>	
Heat-Up Time — Minutes to 80°C	70
Heat-Up Time — Minutes to 150°C	120
Heat-Up Time — Minutes to 220°C	175
Temperature Uniformity	+/- 6% of Set Point
Temperature Stability	+/- 0.2°C @ 150°C
Vacuum Display Range	1100 to 0.1 mTorr
Safety Certifications	UL 61010-1, EN 61010-1, CSA 61010-1, CE



MADE IN  
USA



C US CE






*We reserve the right to alter specifications at any time*

0740590 6/22

Contact Cascade TEK

sales@casca detek.com • casca detek.com • +1-971-371-4096



		<b>Cascade TEK Solutions, LLC</b> 4001 East Plano Parkway, Ste 200 Plano, TX 75074 <a href="http://www.cascadetek.com">www.cascadetek.com</a>	
		NAME: GEB DATE: 10/20/2017	MODEL: <b>TVO-2-SV</b> SIZE: <b>B</b> DWG. NO.: 9100947
THE INFORMATION CONTAINED HEREIN IS THE PROPERTY OF CASCADE TEK SOLUTIONS, LLC. ANY REPRODUCTION OR TRANSMISSION OF THIS INFORMATION IN ANY FORM OR BY ANY MEANS WITHOUT THE WRITTEN PERMISSION OF CASCADE TEK SOLUTIONS, LLC IS PROHIBITED.		SCALE: 1:8 WEIGHT:	REV:

A  
B  
C  
D

A  
B  
C  
D

8 7 6 5 4 3 2 1