

Vacuum Oven

TVO-2-2-SV

PART ID: CTVS223-EA






Capacity	1.67 cu.ft. (47.2 liters)
Interior Dimensions (W x D x H)	12" x 20" x 12" 304 x 508 x 304 mm
Exterior Dimensions (W x D x H)	24.5" x 34.9" x 27.1" 623 x 886 x 689 mm
Unit Weight	183 lb (83 kg)
Shipping Weight	224 lb (102 kg)
Shelves	3 standard, 3 total
Maximum Weight Per Shelf	35 lb / 15.8 kg
Permitted Total Load	105 lb / 47.6 kg
Standard Door Gasket	Silicone Standard
Access Port	1" vacuum port (KF-25)
Controls	Watlow PM9 Plus
Warranty	2 Year Limited Warranty
Interior Construction	300 series stainless steel
Exterior Construction	Painted Cold Roll Steel
Wall Clearance, Each Side	6.0" (152 mm)
Wall Clearance, Rear	12.0" (305 mm)
Vertical Clearance	12.0" (305 mm)
Electrical, 50/60Hz, AC, 1Ø	
Voltage	220V-240V
Full Load Amps	5.5
Nominal Power (watts @150°C)	750
Plug Supplied	EU 1-16P; NEMA 5-15P
Energy Consumption (kWh/day)	
At 80°C	13.9
At 150°C	16.7
At 220°C	26.5
Temperature Range 10°C above ambient to 220°C	
Heat-Up Time — Minutes to 80°C	70
Heat-Up Time — Minutes to 150°C	120
Heat-Up Time — Minutes to 220°C	175
Temperature Uniformity	+/- 6% of Set Point
Temperature Stability	+/- 0.2°C @ 150°C
Vacuum Display Range	1100 to 0.1 mTorr
Safety Certifications	UL 61010-1, EN 61010-1, CSA 61010-1, CE



We reserve the right to alter specifications at any time

0740591 5/24



CASCADE TEK

Cascade TEK Solutions, LLC
300 N 26th Ave Cornelius, OR 97113
www.cascadetek.com

DATE	10/20/2017	MODEL:	TVO-2-2-SV
NAME	GRS	SIZE	DWG. NO. 9100948
DRAWN	GRS	SCALE:	1:8
<small>THE INFORMATION CONTAINED IN THIS SHEET IS UNCLASSIFIED AND IS IN THE PUBLIC DOMAIN. ANY REPRODUCTION OR TRANSMISSION OF THIS INFORMATION IN ANY FORM OR BY ANY MEANS IS PROHIBITED.</small>		REV	

SHEET 1 OF 1