

Vacuum Oven

TVO-5-VC
PART ID: CTVV522



Capacity	4.5 cu.ft. (127.4 liters)
Interior Dimensions (W x D x H)	18.1" x 24.1" x 18.1" 459 x 612 x 459 mm
Exterior Dimensions (W x D x H)	26.6" x 38.7" x 33.2" 676 x 983 x 844 mm
Unit Weight	365 lb (166 kg)
Shipping Weight	460 lb (209 kg)
Shelves	3 standard, 6 total
Maximum Weight Per Shelf	35 lb / 15.8 kg
Permitted Total Load	105 lb / 47.6 kg
Standard Door Gasket	Silicone Standard
Access Port	1" vacuum port (KF-25)
Controls	Watlow F4T
Warranty	2 Year Limited Warranty
Interior Construction	300 series stainless steel
Exterior Construction	Painted Cold Roll Steel
Wall Clearance, Each Side	6.0" (152 mm)
Wall Clearance, Rear	12.0" (305 mm)
Vertical Clearance	12.0" (305 mm)
Electrical, 50/60Hz, AC, 1Ø	
Voltage	110V-120V
Full Load Amps	13
Nominal Power (watts @150°C)	835
Plug Supplied	NEMA 5-15P
Energy Consumption (kWh/day)	
At 80°C	13.2
At 150°C	19.4
At 220°C	30.8
Temperature Range 10°C above ambient to 220°C	
Heat-Up Time — Minutes to 80°C	70
Heat-Up Time — Minutes to 150°C	120
Heat-Up Time — Minutes to 220°C	175
Temperature Uniformity	+/- 6% of Set Point
Temperature Stability	+/- 0.25°C @ 150°C
Operational Vacuum Range	683 to less than 20 mTorr @ 150°C
Vacuum Display Range	1100 to 0.1 mTorr



MADE IN
USA




CE c US



We reserve the right to alter specifications at any time.

0740596 2/20

Contact Cascade TEK
sales@cascadetek.com • cascadetek.com • +1-971-371-4096





Top view


Internal material:
300 series stainless steel

Front view (door not shown)

External material:
Cold rolled mild steel
Powder coated




Anodized Aluminum
Shelf Dimensions



DETAIL B
Shelf standard spacing

Vent port:
1/4" O.D. stainless steel tube



Back view

		Cascade TEK Solutions, LLC 300 N 26th Ave Cornelius, OR 97113 www.cascadetek.com	
		PROPRIETARY AND CONFIDENTIAL THE INFORMATION CONTAINED IN THIS DRAWING IS THE SOLE PROPERTY OF CASCADE TEK SOLUTIONS, LLC. IT IS TO BE USED FOR THE PART OR AS A WHOLE WITHOUT THE WRITTEN PERMISSION OF Cascade TEK SOLUTIONS, LLC.	NAME: GFB DATE: 10/24/2017
MODEL: TWO-5-VC		SIZE: B DWG. NO.: 9100943	REV:
SCALE: 1:11		WEIGHT:	SHEET 1 OF 1