


Vacuum Oven

TVO-5-2-VC
PART ID: CTVV522-EA



Capacity	4.5 cu.ft. (127.4 liters)
Interior Dimensions (W x D x H)	18.1" x 24.1" x 18.1" 459 x 612 x 459 mm
Exterior Dimensions (W x D x H)	26.6" x 38.7" x 33.2" 676 x 983 x 844 mm
Unit Weight	365 lb (166 kg)
Shipping Weight	460 lb (209 kg)
Shelves	3 standard, 6 total
Maximum Weight Per Shelf	35 lb / 15.8 kg
Permitted Total Load	105 lb / 47.6 kg
Standard Door Gasket	Silicone Standard
Access Port	1" vacuum port (KF-25)
Controls	Watlow F4T
Warranty	2 Year Limited Warranty
Interior Construction	300 series stainless steel
Exterior Construction	Painted Cold Roll Steel
Wall Clearance, Each Side	6.0" (152 mm)
Wall Clearance, Rear	12.0" (305 mm)
Vertical Clearance	12.0" (305 mm)
Electrical, 50/60Hz, AC, 1Ø	
Voltage	220V-240V
Full Load Amps	7.0
Nominal Power (watts @150°C)	936
Plug Supplied	EU 1-16P; NEMA 6-15P
Energy Consumption (kWh/day)	
At 80°C	18.1
At 150°C	21.8
At 220°C	34.5
Temperature Range 10°C above ambient to 220°C	
Heat-Up Time — Minutes to 80°C	70
Heat-Up Time — Minutes to 150°C	120
Heat-Up Time — Minutes to 220°C	175
Temperature Uniformity	+/- 6% of Set Point
Temperature Stability	+/- 0.25°C @ 150°C
Operational Vacuum Range	683 to less than 20 mTorr @ 150°C
Vacuum Display Range	1100 to 0.1 mTorr



MADE IN
USA





CE C US



We reserve the right to alter specifications at any time





0740597 2/20



External material:
Cold rolled mild steel
Powder coated

Internal material:
300 series stainless steel

External material:
Cold rolled mild steel
Powder coated



		Cascade TEK Solutions, LLC 300 N 26th Ave Cornelius, OR 97113 www.cascadetek.com	
PROPRIETARY AND CONFIDENTIAL	DATE: 10/26/2017	MODEL: T20-5-2-VC	REV: 1
NAME: GBS	DATE: 10/26/2017	SIZE: B	DWG. NO.: 9100958
THE INFORMATION CONTAINED IN THIS DRAWING IS THE PROPERTY OF CASCADE TEK SOLUTIONS, LLC. ANY REPRODUCTION OR TRANSMISSION IN ANY FORM OR BY ANY MEANS WITHOUT THE WRITTEN PERMISSION OF CASCADE TEK SOLUTIONS, LLC IS PROHIBITED.		SCALE: 1:1	WEIGHT: SHEET 1 OF 1