

Vacuum Oven


SVO-5-VC

PART ID: CTVR522-H



Capacity	4.88 cu.ft. (138 liters)
Interior Dimensions (W x D x H)	18" x 26" x 18" 457 x 660 x 457 mm
Exterior Dimensions (W x D x H)	28.8" x 46.5" x 64.6" 732 x 1182 x 1640 mm
Unit Weight	760 lb (345 kg)
Shipping Weight	1078 lb (489 kg)
Shelves	3 Anodized Aluminum Standard, 6 Maximum
Maximum Weight Per Shelf	50 lb / 22.7 kg
Permitted Total Load	200 lb / 90.7 kg
Standard Door Gasket	Viton Standard
Access Port	1-KF25
Controls	Watlow F4T
Warranty	2 Year Limited Warranty
Interior Construction	300 series stainless steel
Exterior Construction	20 gauge steel, powder coated
Wall Clearance, Each Side	6.0" (152 mm)
Wall Clearance, Rear	12.0" (305 mm)
Vertical Clearance	12.0" (305 mm)
Electrical, 50/60Hz, AC, 1Ø	
Voltage	220V-240V
Power Current (Max Draw) Amps	20
Nominal Power (watts @150°C)	680
Direct Hardwired	✓
Energy Consumption (kWh/day)	
At 80°C	9.22
At 150°C	16.34
At 220°C	26.3
Temperature Range	
	Ambient 10°C to 220°C
Heat-Up Time — Minutes to 80°C	68
Heat-Up Time — Minutes to 150°C	125
Heat-Up Time — Minutes to 220°C	177
Operation Vacuum Range	
	760 Torr to Less than 10 mTorr
Temperature Uniformity	+/- 6% of Set Point
Temperature Stability	+/- 0.2°C @ 150°C
Vacuum Display Range	1100 to 0.1 mTorr


We reserve the right to alter specifications at any time.




0740648 5/21

Contact Cascade TEK


sales@cascadetek.com • cascadetek.com • +1-971-371-4096




OPTIONAL CASTERS ADD
0.75 IN TO OVERALL HEIGHT
DETAIL A
SCALE 1 : 4




Top View




EXTERIOR MATERIAL:
COLD ROLLED STEEL
POWDER COATED



Back View



Side View



INTERNAL MATERIAL:
300 SERIES STAINLESS STEEL

Front View (Door Not Shown)

CASCADE TEK				Cascade TEK Solutions, Inc. 2010 Century Center Blvd, Ste 21 Irving, Texas 75062 www.cascadetek.com			
PROPRIETARY AND CONFIDENTIAL THE INFORMATION CONTAINED IN THIS DRAWING IS THE SOLE PROPERTY OF CASCADE TEK SOLUTIONS, INC. AND IS TO BE USED ONLY FOR THE PROJECT AND PART OR AS A WHOLE WITHOUT THE WRITTEN PERMISSION OF Cascade TEK SOLUTIONS, INC.				MODEL: SVO-5-VC			
DRAWN	NAME	DATE	SCALE	SIZE	DWG. NO.	REV	
GT	GT	1/15/2020	SCALE: 1:20	B	9101005	A	
				WEIGHT: SHEET 1 OF 4			